



# ENOX1 Antibody, HRP conjugated

Product Code	CSB-PA007676LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TC92
Immunogen	Recombinant Human Ecto-NOX disulfide-thiol exchanger 1 protein (101-340AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ecto-NOX disulfide-thiol exchanger 1 (Candidate growth-related and time keeping constitutive hydroquinone [NADH] oxidase) (cCNOX) (Cell proliferation-inducing gene 38 protein) (Constitutive Ecto-NOX) (cNOX) [Includes: Hydroquinone [NADH] oxidase (EC 1.-.-.-); Protein disulfide-thiol oxidoreductase (EC 1.-.-.-)], ENOX1
Immunogen Species	Homo sapiens (Human)
Research Area	Cancer
Target Names	ENOX1