



EPO Antibody

Product Code	CSB-PA007743LA01MOV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P07865
Immunogen	Recombinant Macaca fascicularis Erythropoietin protein (28-192AA)
Raised In	Rabbit
Species Reactivity	Macaca fascicularis
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Erythropoietin, EPO
Immunogen Species	Macaca fascicularis (Crab-eating macaque) (Cynomolgus monkey)
Research Area	Others
Target Names	EPO
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.