



EPO Antibody, FITC conjugated

Product Code	CSB-PA007743LC01MOV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P07865
Immunogen	Recombinant Macaca fascicularis Erythropoietin protein (28-192AA)
Raised In	Rabbit
Species Reactivity	Macaca fascicularis
Tested Applications	ELISA
Relevance	Erythropoietin is the principal hormone involved in the regulation of erythrocyte differentiation and the maintenance of a physiological level of circulating erythrocyte mass.
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Erythropoietin, EPO
Species	Macaca fascicularis (Crab-eating macaque) (Cynomolgus monkey)
Research Area	Others
Target Names	EPO