



ERG Antibody, FITC conjugated

Product Code	CSB-PA007781HC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P11308
Immunogen	Recombinant Human Transcriptional regulator ERG protein (1-486AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transcriptional regulator ERG (Transforming protein ERG), ERG
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	ERG
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.