



F11 Antibody, HRP conjugated

Product Code	CSB-PA007916LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P03951
Immunogen	Recombinant Human Coagulation factor XI protein (19-387AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Coagulation factor XI (FXI) (EC 3.4.21.27) (Plasma thromboplastin antecedent) (PTA) [Cleaved into: Coagulation factor XIa heavy chain; Coagulation factor XIa light chain], F11
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	F11
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.