




FAH Antibody

Product Code	CSB-PA007965LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P16930
Immunogen	Recombinant Human Fumarylacetoacetase protein (2-250AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC, IF; Recommended dilution: WB:1:2000-1:5000, IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Fumarylacetoacetase (FAA) (EC 3.7.1.2) (Beta-diketonase) (Fumarylacetoacetate hydrolase), FAH
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	FAH

Image



Western Blot

Positive WB detected in: Mouse kidney tissue, Mouse stomach tissue


All lanes: FAH antibody at 3.4µg/ml

Secondary


Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 47, 39 kDa

Observed band size: 47 kDa



Immunohistochemistry of paraffin-embedded human liver tissue using CSB-PA007965LA01HU at dilution of 1:100



Immunofluorescent analysis of HeLa cells using CSB-PA007965LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)