



FAM107A Antibody, FITC conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA007987LC01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | O95990 |
| Immunogen | Recombinant Human Actin-associated protein FAM107A protein (1-144AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Form | Liquid |
| Conjugate | FITC |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Actin-associated protein FAM107A (Down-regulated in renal cell carcinoma 1) (Protein TU3A), FAM107A, DRR1 TU3A |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Cell Biology |
| Target Names | FAM107A |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |