



NUTM2F Antibody

Product Code	CSB-PA008184LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A1L443
Immunogen	Recombinant Human NUT family member 2F protein (450-690AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	NUT family member 2F, NUTM2F, FAM22F
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	NUTM2F
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.