



FFAR2 Antibody, HRP conjugated

Product Code	CSB-PA008605LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15552
Immunogen	Recombinant Human Free fatty acid receptor 2 protein (277-330AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Free fatty acid receptor 2 (G-protein coupled receptor 43), FFAR2, FFA2 GPCR43 GPR43
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	FFAR2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.