



FFAR2 Antibody, FITC conjugated

Product Code	CSB-PA008605LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15552
Immunogen	Recombinant Human Free fatty acid receptor 2 protein (277-330AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Free fatty acid receptor 2 (G-protein coupled receptor 43), FFAR2, FFA2 GPCR43 GPR43
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	FFAR2