



FFAR3 Antibody, FITC conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA008606LC01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | O14843 |
| Immunogen | Recombinant Human Free fatty acid receptor 3 protein (280-346AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Form | Liquid |
| Conjugate | FITC |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Free fatty acid receptor 3 (G-protein coupled receptor 41), FFAR3, GPR41 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Cardiovascular |
| Target Names | FFAR3 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |