



FFAR3 Antibody, Biotin conjugated

Product Code	CSB-PA008606LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O14843
Immunogen	Recombinant Human Free fatty acid receptor 3 protein (280-346AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Free fatty acid receptor 3 (G-protein coupled receptor 41), FFAR3, GPR41
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	FFAR3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.