




G6PC Antibody

Product Code	CSB-PA009118LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P35575
Immunogen	Recombinant Human Glucose-6-phosphatase protein (82-117AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glucose-6-phosphatase (G-6-Pase) (G6Pase) (EC 3.1.3.9) (Glucose-6-phosphatase alpha) (G6Pase-alpha), G6PC, G6PT
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	G6PC

Image



Immunohistochemistry of paraffin-embedded human liver tissue using CSB-PA009118LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.