





GGTLC1 Antibody

Product Code	CSB-PA009404LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BX51
Immunogen	Recombinant Human Glutathione hydrolase light chain 1 protein (1-225AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glutathione hydrolase light chain 1 (Gamma-glutamyltransferase light chain 1) (Gamma-glutamyltransferase-like activity 4) (Gamma-glutamyltransferase-like protein 6), GGTLC1, GGTLA4
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	GGTLC1

Image



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA009404LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human pancreatic cancer using CSB-PA009404LA01HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.