



GLS Antibody, Biotin conjugated

Product Code	CSB-PA009528LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O94925
Immunogen	Recombinant Human Glutaminase kidney isoform, mitochondrial protein (616-669AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glutaminase kidney isoform, mitochondrial (GLS) (EC 3.5.1.2) (K-glutaminase) (L-glutamine amidohydrolase), GLS, GLS1 KIAA0838
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	GLS
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.