



# GLTP Antibody

Product Code	CSB-PA009539LA01BO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P68265
Immunogen	Recombinant Bovine Glycolipid transfer protein (2-209AA)
Raised In	Rabbit
Species Reactivity	Bovine
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glycolipid transfer protein (GLTP), GLTP
Immunogen Species	Bos taurus (Bovine)
Research Area	Others
Target Names	GLTP