



GLTP Antibody

Product Code	CSB-PA009539LA01BO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P68265
Immunogen	Recombinant Bovine Glycolipid transfer protein (2-209AA)
Raised In	Rabbit
Species Reactivity	Bovine
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glycolipid transfer protein (GLTP), GLTP
Immunogen Species	Bos taurus (Bovine)
Research Area	Others
Target Names	GLTP
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.