



GLTP Antibody, FITC conjugated

Product Code	CSB-PA009539LC01BO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P68265
Immunogen	Recombinant Bovine Glycolipid transfer protein (2-209AA)
Raised In	Rabbit
Species Reactivity	Bovine
Tested Applications	ELISA
Relevance	Accelerates the intermembrane transfer of various glycolipids. Catalyzes the transfer of various glycosphingolipids between membranes but does not catalyze the transfer of phospholipids. May be involved in the intracellular translocation of glucosylceramides.
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glycolipid transfer protein (GLTP), GLTP
Immunogen Species	Bos taurus (Bovine)
Research Area	Others
Target Names	GLTP