



GMFG Antibody

Product Code	CSB-PA009574LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60234
Immunogen	Recombinant Human Glia maturation factor gamma protein (1-142AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glia maturation factor gamma (GMF-gamma), GMFG
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	GMFG
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.