



GNG11 Antibody, Biotin conjugated

Product Code	CSB-PA009611LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P61952
Immunogen	Recombinant Human Guanine nucleotide-binding protein G(I)/G(S)/G(O) subunit gamma-11 protein (1-70AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Guanine nucleotide-binding protein G(I)/G(S)/G(O) subunit gamma-11, GNG11, GNGT11
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	GNG11
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.