



GRB2 Antibody, Biotin conjugated

Product Code	CSB-PA009889ED01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P62993
Immunogen	Recombinant Human Growth factor receptor-bound protein 2 protein (11-200AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Growth factor receptor-bound protein 2 (Adapter protein GRB2) (Protein Ash) (SH2/SH3 adapter GRB2), GRB2, ASH
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	GRB2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.