



GSK3B Antibody, HRP conjugated

Product Code	CSB-PA009963LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P49841
Immunogen	Recombinant Human Glycogen synthase kinase-3 beta protein (1-420AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glycogen synthase kinase-3 beta (GSK-3 beta) (EC 2.7.11.26) (Serine/threonine-protein kinase GSK3B) (EC 2.7.11.1), GSK3B
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	GSK3B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.