



GUCA1A Antibody

Product Code	CSB-PA010044LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P43080
Immunogen	Recombinant Human Guanylyl cyclase-activating protein 1 protein (2-201AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Guanylyl cyclase-activating protein 1 (GCAP 1) (Guanylate cyclase activator 1A), GUCA1A, C6orf131 GCAP GCAP1 GUCA1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	GUCA1A