



GUCA1A Antibody

Product Code	CSB-PA010044LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P43080
Immunogen	Recombinant Human Guanylyl cyclase-activating protein 1 protein (2-201AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Guanylyl cyclase-activating protein 1 (GCAP 1) (Guanylate cyclase activator 1A), GUCA1A, C6orf131 GCAP GCAP1 GUCA1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	GUCA1A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.