



HEXIM2 Antibody, HRP conjugated

Product Code	CSB-PA010319LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96MH2
Immunogen	Recombinant Human Protein HEXIM2 protein (1-286AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein HEXIM2 (Hexamethylene bis-acetamide-inducible protein 2), HEXIM2
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	HEXIM2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.