




HS3ST1 Antibody

Product Code	CSB-PA010747LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O14792
Immunogen	Recombinant Human Heparan sulfate glucosamine 3-O-sulfotransferase 1 protein (21-307AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Heparan sulfate glucosamine 3-O-sulfotransferase 1 (EC 2.8.2.23) (Heparan sulfate D-glucosaminyl 3-O-sulfotransferase 1) (3-OST-1) (Heparan sulfate 3-O-sulfotransferase 1) (h3-OST-1), HS3ST1, 3OST 3OST1
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	HS3ST1

Image



Immunohistochemistry of paraffin-embedded human heart tissue using CSB-PA010747LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.