



IFI44 Antibody, FITC conjugated

Product Code	CSB-PA011014LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TCB0
Immunogen	Recombinant Human Interferon-induced protein 44 protein (101-300AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Interferon-induced protein 44 (p44) (Microtubule-associated protein 44), IFI44, MTAP44
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	IFI44
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.