



IFITM2 Antibody, Biotin conjugated

Product Code	CSB-PA011025LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q01629
Immunogen	Recombinant Human Interferon-induced transmembrane protein 2 protein (1-56AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Interferon-induced transmembrane protein 2 (Dispanin subfamily A member 2c) (DSPA2c) (Interferon-inducible protein 1-8D), IFITM2
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	IFITM2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.