



Irf5 Antibody

| | |
|----------------------------|--|
| Product Code | CSB-PA011820LA01MO |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P56477 |
| Immunogen | Recombinant Mouse Interferon regulatory factor 5 protein (1-497AA) |
| Raised In | Rabbit |
| Species Reactivity | Mouse |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Interferon regulatory factor 5 (IRF-5), Irf5 |
| Immunogen Species | Mus musculus (Mouse) |
| Research Area | Immunology |
| Target Names | Irf5 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |