



ITGAL Antibody, HRP conjugated

Product Code	CSB-PA011875LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P20701
Immunogen	Recombinant Human Integrin alpha-L protein (201-500AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Integrin alpha-L (CD11 antigen-like family member A) (Leukocyte adhesion glycoprotein LFA-1 alpha chain) (LFA-1A) (Leukocyte function-associated molecule 1 alpha chain) (CD antigen CD11a), ITGAL, CD11A
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ITGAL
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.