



ITM2A Antibody, HRP conjugated

Product Code	CSB-PA011903LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43736
Immunogen	Recombinant Human Integral membrane protein 2A protein (79-263AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Integral membrane protein 2A (Protein E25), ITM2A
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	ITM2A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.