




KIAA0513 Antibody

Product Code	CSB-PA012195LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60268
Immunogen	Recombinant Human Uncharacterized protein KIAA0513 protein (1-301AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IF; Recommended dilution: WB:1:2000-1:5000, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Uncharacterized protein KIAA0513, KIAA0513
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	KIAA0513

Image



Western Blot

Positive WB detected in: Jurkat whole cell lysate, K562 whole cell lysate, Mouse brain tissue


All lanes: KIAA0513 antibody at 2.7µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 47, 34, 46 kDa

Observed band size: 57 kDa



Immunofluorescent analysis of HepG2 cells using CSB-PA012195LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.