




KIAA1324L Antibody

Product Code	CSB-PA012251LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A8MWY0
Immunogen	Recombinant Human UPF0577 protein KIAA1324-like protein (227-526AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	UPF0577 protein KIAA1324-like (Estrogen-induced gene 121-like protein) (hEIG121L), KIAA1324L, EIG121L
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	KIAA1324L

Image



Immunohistochemistry of paraffin-embedded human gastric cancer using CSB-PA012251LA01HU at dilution of 1:100

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.