



KLC4 Antibody

Product Code	CSB-PA012382LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NSK0
Immunogen	Recombinant Human Kinesin light chain 4 protein (401-619AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Kinesin light chain 4 (KLC 4) (Kinesin-like protein 8), KLC4, KNSSL8
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	KLC4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.