



# Klkb1 Antibody, FITC conjugated

<b>Product Code</b>	CSB-PA012461YC01RA
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P14272
<b>Immunogen</b>	Recombinant Rat Plasma kallikrein protein (391-638AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Rat
<b>Form</b>	Liquid
<b>Conjugate</b>	FITC
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Plasma kallikrein (EC 3.4.21.34) (Fletcher factor) (Kininogenin) (Plasma prekallikrein) [Cleaved into: Plasma kallikrein heavy chain; Plasma kallikrein light chain], Klkb1, Klk3 Pk
<b>Immunogen Species</b>	Rattus norvegicus (Rat)
<b>Research Area</b>	Cardiovascular
<b>Target Names</b>	Klkb1
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.