



Klkb1 Antibody, Biotin conjugated

Product Code	CSB-PA012461YD01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P14272
Immunogen	Recombinant Rat Plasma kallikrein protein (391-638AA)
Raised In	Rabbit
Species Reactivity	Rat
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Plasma kallikrein (EC 3.4.21.34) (Fletcher factor) (Kininogenin) (Plasma prekallikrein) [Cleaved into: Plasma kallikrein heavy chain; Plasma kallikrein light chain], Klkb1, Klk3 Pk
Immunogen Species	Rattus norvegicus (Rat)
Research Area	Cardiovascular
Target Names	Klkb1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.