



KRT17 Antibody, Biotin conjugated

Product Code	CSB-PA012516LD01BO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A1L595
Immunogen	Recombinant Bovine Keratin, type I cytoskeletal 17 protein (1-441AA)
Raised In	Rabbit
Species Reactivity	Bovine
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Keratin, type I cytoskeletal 17 (Cytokeratin-17) (CK-17) (Keratin-17) (K17), KRT17
Immunogen Species	Bos taurus (Bovine)
Research Area	Others
Target Names	KRT17
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.