



KRT81 Antibody, Biotin conjugated

Product Code	CSB-PA012576LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14533
Immunogen	Recombinant Human Keratin, type II cuticular Hb1 protein (3-505AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Keratin, type II cuticular Hb1 (Hair keratin K2.9) (Keratin, hair, basic, 1) (Keratin-81) (K81) (Metastatic lymph node 137 gene protein) (MLN 137) (Type II hair keratin Hb1) (Type-II keratin Kb21) (ghHKb1) (ghHb1), KRT81, KRTHB1 MLN137
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	KRT81
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.