



L2HGDH Antibody, Biotin conjugated

Product Code	CSB-PA012709LD01BO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A7MBI3
Immunogen	Recombinant Bovine L-2-hydroxyglutarate dehydrogenase, mitochondrial protein (53-463AA)
Raised In	Rabbit
Species Reactivity	Bovine
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	L-2-hydroxyglutarate dehydrogenase, mitochondrial (EC 1.1.99.2) (Duranin), L2HGDH
Immunogen Species	Bos taurus (Bovine)
Research Area	Signal Transduction
Target Names	L2HGDH
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.