



LHFPL6 Antibody, Biotin conjugated

Product Code	CSB-PA012912LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y693
Immunogen	Recombinant Human LHFPL tetraspan subfamily member 6 protein (31-83AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	LHFPL tetraspan subfamily member 6 protein (Lipoma HMGIC fusion partner), LHFPL6, LHFP
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	LHFPL6