



LHFPL6 Antibody, Biotin conjugated

Product Code	CSB-PA012912LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y693
Immunogen	Recombinant Human LHFPL tetraspan subfamily member 6 protein (31-83AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	LHFPL tetraspan subfamily member 6 protein (Lipoma HMGIC fusion partner), LHFPL6, LHFP
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	LHFPL6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.