



LIMS3 Antibody, Biotin conjugated

Product Code	CSB-PA012957LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0CW19
Immunogen	Recombinant Human LIM and senescent cell antigen-like-containing domain protein 3 protein (1-117AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	LIM and senescent cell antigen-like-containing domain protein 3 (Particularly interesting new Cys-His protein 3) (PINCH-3), LIMS3, PINCH3
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	LIMS3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.