



LIN7A Antibody, HRP conjugated

Product Code	CSB-PA012963LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O14910
Immunogen	Recombinant Human Protein lin-7 homolog A protein (1-233AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein lin-7 homolog A (Lin-7A) (hLin-7) (Mammalian lin-seven protein 1) (MALS-1) (Tax interaction protein 33) (TIP-33) (Vertebrate lin-7 homolog 1) (Veli-1), LIN7A, MALS1 VELI1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	LIN7A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.