



LMF2 Antibody, HRP conjugated

Product Code	CSB-PA013001LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BU23
Immunogen	Recombinant Human Lipase maturation factor 2 protein (432-636AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Lipase maturation factor 2 (Transmembrane protein 112B) (Transmembrane protein 153), LMF2, TMEM112B TMEM153
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	LMF2