



LMF2 Antibody, Biotin conjugated

Product Code	CSB-PA013001LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BU23
Immunogen	Recombinant Human Lipase maturation factor 2 protein (432-636AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Lipase maturation factor 2 (Transmembrane protein 112B) (Transmembrane protein 153), LMF2, TMEM112B TMEM153
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	LMF2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.