



LPAR4 Antibody, Biotin conjugated

Product Code	CSB-PA013050LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99677
Immunogen	Recombinant Human Lysophosphatidic acid receptor 4 protein (1-43AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Lysophosphatidic acid receptor 4 (LPA receptor 4) (LPA-4) (G-protein coupled receptor 23) (P2Y purinoceptor 9) (P2Y9) (P2Y5-like receptor) (Purinergic receptor 9), LPAR4, GPR23 LPA4 P2RY9
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	LPAR4