



LPAR4 Antibody, Biotin conjugated

Product Code	CSB-PA013050LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99677
Immunogen	Recombinant Human Lysophosphatidic acid receptor 4 protein (1-43AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Lysophosphatidic acid receptor 4 (LPA receptor 4) (LPA-4) (G-protein coupled receptor 23) (P2Y purinoceptor 9) (P2Y9) (P2Y5-like receptor) (Purinergic receptor 9), LPAR4, GPR23 LPA4 P2RY9
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	LPAR4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.