



LRRN2 Antibody, HRP conjugated

Product Code	CSB-PA013183LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75325
Immunogen	Recombinant Human Leucine-rich repeat neuronal protein 2 protein (201-500AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Leucine-rich repeat neuronal protein 2 (Glioma amplified on chromosome 1 protein) (Leucine-rich repeat neuronal protein 5), LRRN2, GAC1 LRRN5
Immunogen Species	Homo sapiens (Human)
Research Area	Cancer
Target Names	LRRN2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.