



LY6G6C Antibody

Product Code	CSB-PA013245LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95867
Immunogen	Recombinant Human Lymphocyte antigen 6 complex locus protein G6c protein (19-99AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Lymphocyte antigen 6 complex locus protein G6c, LY6G6C, C6orf24 G6C NG24
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	LY6G6C
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.