



LYPD3 Antibody, Biotin conjugated

Product Code	CSB-PA013263LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95274
Immunogen	Recombinant Human Ly6/PLAUR domain-containing protein 3 protein (31-326AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ly6/PLAUR domain-containing protein 3 (GPI-anchored metastasis-associated protein C4.4A homolog) (Matrigel-induced gene C4 protein) (MIG-C4), LYPD3, C4.4A
Immunogen Species	Homo sapiens (Human)
Research Area	Cancer
Target Names	LYPD3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.